Packet Types OverView

Type ID	Type Name	From	То	Type Description
0	Scene Packet	brachlOplexus	Unity	Contains the scene index, frame rate, and initial camera positions and joint limits for scene initialization in Unity
1	Control Packet	brachlOplexus	Unity	Contains the ID, motor state and velocity of each joint in brachIOplexus
2	<u>Update Packet</u>	brachlOplexus	Unity	Contains the updated camera positions and joint limits for a running scene
3	Camera Request Packet	brachlOplexus	Unity	An empty packet that acts as a request to retrieve current camera positions from Unity
4	Feedback Packet	Unity	brachlOplexus	Contains the ID, current position and current velocity of each joint in Unity
5	Camera Position Packet	Unity	brachlOplexus	Contains the current camera positions in Unity

Scene Packet (Type 0) Structure

Scene Packet (Type 0) Structure		
Byte Number	Info	Value
0	HEADER	255
1	HEADER	255
2	Packet Type ID	0
3	Length of Data	40
4	Scene Index	1 to the number of total scenes
5	Frame Rate	Any positive integer
6	Camera Position X - low byte	
7	Camera Position X - high byte	
8	Camera Position X - sign indicator	1 if position x is positive, 0 if position x is negative
9	Camera Position Y - low byte	
10	Camera Position Y - high byte	
11	Camera Position Y - sign indicator	1 if position y is positive, 0 if position y is negative
12	Camera Position Z - low byte	
13	Camera Position Z - high byte	
14	Camera Position Z - sign indicator	1 if position z is positive, 0 if position z is negative
15	Camera Rotation X - low byte	
16	Camera Rotation X - high byte	
17	Camera Rotation X - sign indicator	1 if rotation x is positive, 0 if rotation x is negative
18	Camera Rotation Y - low byte	
19	Camera Rotation Y - high byte	
20	Camera Rotation Y - sign indicator	1 if rotation y is positive, 0 if rotation y is negative
21	Camera Rotation Z - low byte	
22	Camera Rotation Z - high byte	
23	Camera Rotation Z - sign indicator	1 if rotation z is positive, 0 if rotation z is negative
24	Shoulder Minimum Position - low byte	
25	Shoulder Minimum Position - high byte	
26	Shoulder Maximum Position - low byte	
27	Shoulder Maximum Position - high byte	
28	Elbow Minimum Position - low byte	
29	Elbow Minimum Position - high byte	
30	Elbow Maximum Position - low byte	
31	Elbow Maximum Position - high byte	
32	Wrist Rotate Minimum Position - low byte	
33	Wrist Rotate Minimum Position - high byte	
34	Wrist Rotate Maximum Position - low byte	
35	Wrist Rotate Maximum Position - high byte	
36	Wrist Flex Minimum Position - low byte	
37	Wrist Flex Minimum Position - high byte	
38	Wrist Flex Maximum Position - low byte	
39	Wrist Flex Maximum Position - high byte	
40	Hand Minimum Position - low byte	
41	Hand Minimum Position - high byte	
42	Hand Maximum Position - low byte	
43	Hand Maximum Position - high byte	
44	Checksum	

Control Packet (Type 1) Structure

Byte Number	Info	Value
0	HEADER	255
1	HEADER	255
2	Packet Type ID	1
3	Length of Data	20
4	Shoulder Joint ID	0
5	Shoulder Joint Velocity - low byte	
6	Shoulder Joint Velocity - high byte	
7	Shoulder Joint Motor State	0~3
8	Elbow Joint ID	1
9	Elbow Joint Velocity - low byte	
10	Elbow Joint Velocity - high byte	
11	Elbow Joint Motor State	0~3
12	Wrist Rotate Joint ID	2
13	Wrist Rotate Joint Velocity - low byte	
14	Wrist Rotate Joint Velocity - high byte	
15	Wrist Rotate Joint Motor State	0~3
16	Wrist Flex Joint ID	3
17	Wrist Flex Joint Velocity - low byte	
18	Wrist Flex Joint Velocity - high byte	
19	Wrist Flex Joint Motor State	0~3
20	Hand Joint ID	4
21	Hand Joint Velocity - low byte	
22	Hand Joint Velocity - high byte	
23	Hand Joint Motor State	0~3
24	Checksum	

Update Packet (Type 2) Structure

Update Packet (Type 2) Structure				
Byte Number	Info	Value		
0	HEADER	255		
1	HEADER	255		
2	Packet Type ID	2		
3	Length of Data	38		
4	Camera Position X - low byte			
5	Camera Position X - high byte			
6	Camera Position X - sign indicator	1 if position x is positive, 0 if position x is negative		
7	Camera Position Y - low byte			
8	Camera Position Y - high byte			
9	Camera Position Y - sign indicator	1 if position y is positive, 0 if position y is negative		
10	Camera Position Z - low byte			
11	Camera Position Z - high byte			
12	Camera Position Z - sign indicator	1 if position z is positive, 0 if position z is negative		
13	Camera Rotation X - low byte			
14	Camera Rotation X - high byte			
15	Camera Rotation X - sign indicator	1 if rotation x is positive, 0 if rotation x is negative		
16	Camera Rotation Y - low byte			
17	Camera Rotation Y - high byte			
18	Camera Rotation Y - sign indicator	1 if rotation y is positive, 0 if rotation y is negative		
19	Camera Rotation Z - low byte			
20	Camera Rotation Z - high byte			
21	Camera Rotation Z - sign indicator	1 if rotation z is positive, 0 if rotation z is negative		
22	Shoulder Minimum Position - low byte			
23	Shoulder Minimum Position - high byte			
24	Shoulder Maximum Position - low byte			
25	Shoulder Maximum Position - high byte			
26	Elbow Minimum Position - low byte			
27	Elbow Minimum Position - high byte			
28	Elbow Maximum Position - low byte			
29	Elbow Maximum Position - high byte			
30	Wrist Rotate Minimum Position - low byte			
31	Wrist Rotate Minimum Position - high byte			
32	Wrist Rotate Maximum Position - low byte			
33	Wrist Rotate Maximum Position - high byte			
34	Wrist Flex Minimum Position - low byte			
35	Wrist Flex Minimum Position - high byte			
36	Wrist Flex Maximum Position - low byte			
37	Wrist Flex Maximum Position - high byte			
38	Hand Minimum Position - low byte			
39	Hand Minimum Position - high byte			
40	Hand Maximum Position - low byte			
41	Hand Maximum Position - high byte			
42	Checksum			

Camera Request Packet (Type 3) Structure

Byte Number	Info	Value
0	HEADER	255
1	HEADER	255
2	Packet Type ID	3
3	Checksum	

Feedback Packet (Type 4) Structure

Feedback Packet (Type 4) Structure				
Byte Number	Info	Value		
0	HEADER	255		
1	HEADER	255		
2	Packet Type ID	4		
3	Length of Data	30		
4	Shoulder Joint ID	0		
5	Shoulder Joint Position - low byte			
6	Shoulder Joint Position - high byte			
7	Shoulder Joint Position - sign indicator	1 if joint position is positive, 0 if joint position is negative		
8	Shoulder Joint Velocity - low byte			
9	Shoulder Joint Velocity - high byte			
10	Elbow Joint ID	1		
11	Elbow Joint Position - low byte			
12	Elbow Joint Position - high byte			
13	Elbow Joint Position - sign indicator	1 if joint position is positive, 0 if joint position is negative		
14	Elbow Joint Velocity - low byte			
15	Elbow Joint Velocity - high byte			
16	Wrist Rotate Joint ID	2		
17	Wrist Rotate Joint Position - low byte			
18	Wrist Rotate Joint Position - high byte			
19	Wrist Rotate Joint Position - sign indicator	1 if joint position is positive, 0 if joint position is negative		
20	Wrist Rotate Joint Velocity - low byte			
21	Wrist Rotate Joint Velocity - high byte			
22	Wrist Flex Joint ID	3		
23	Wrist Flex Joint Position - low byte			
24	Wrist Flex Joint Position - high byte			
25	Wrist Flex Joint Position - sign indicator	1 if joint position is positive, 0 if joint position is negative		
26	Wrist Flex Joint Velocity - low byte			
27	Wrist Flex Joint Velocity - high byte			
28	Hand Joint ID	4		
29 Hand Joint Position - low byte				
30	Hand Joint Position - high byte			
31	Hand Joint Position - sign indicator	1 if joint position is positive, 0 if joint position is negative		
32	Hand Joint Velocity - low byte			
33	Hand Joint Velocity - high byte			
34	Checksum			

Camera Position Packet (Type 5) Structure

Camera Position Packet (Type 5) Structure			
Byte Number	Info	Value	
0	HEADER	255	
1	HEADER	255	
2	Packet Type ID	5	
3	Length of Data	18	
4	Camera Position X - low byte		
5	Camera Position X - high byte		
6	Camera Position X - sign indicator	1 if position x is positive, 0 if position x is negative	
7	Camera Position Y - low byte		
8	Camera Position Y - high byte		
9	Camera Position Y - sign indicator	1 if position y is positive, 0 if position y is negative	
10	Camera Position Z - low byte		
11	Camera Position Z - high byte		
12	Camera Position Z - sign indicator	1 if position z is positive, 0 if position z is negative	
13	Camera Rotation X - low byte		
14	Camera Rotation X - high byte		
15	Camera Rotation X - sign indicator	1 if rotation x is positive, 0 if rotation x is negative	
16	Camera Rotation Y - low byte		
17	Camera Rotation Y - high byte		
18	Camera Rotation Y - sign indicator	1 if rotation y is positive, 0 if rotation y is negative	
19	Camera Rotation Z - low byte		
20	Camera Rotation Z - high byte		
21	Camera Rotation Z - sign indicator	1 if rotation z is positive, 0 if rotation z is negative	
22	Checksum		